



BOMB CALORIMETER



AUTOMATIC BOMB CALORIMETER (DIGITAL) (WITH SAFETY DEVICE)

Our Bomb Calorimeter is widely used for determining gross calorific values of liquid and solid samples, Automatically. User Doesn't need to make any manual calculation, all calculation done automatically by the system itself & then are Automatically Displayed, Printed & Transformed to PC.

All the processor is done by highly intelligent Microprocessor based program. We use special imported calibrators for proper calibration of the system.

The controllers are specially designed & manufactured to achieve best possible measuring accuracy desired by the user. **Supplied without oxygen cylinder**

SPECIAL FEATURES:

- **Fully Automatic Micro Processor based.**
- **Automatic measurement and calculation of CALORIFIC VALUE/WATER EQUIVALENT.**
- **BIG Jumbo LCD Display (BLUE)**
- **In built Printer.**
- **16 soft touch keypad.**
- **Weight of tablet feedable.**
- **Water equivalent feedable.**
- **Full test report printout with data and time & company name.**
- **Temperature Scanning Resolution: 0.01°C.**
- **Alarm on Firing & after the test complete.**
- **Sensor Open detection.**
- **Fuse wire open detection.**
- **No paper or printer door open detection.**
- **User interaction message on LCD.**
- **Type of test Selectable.**
- **Memory storage for data.**
- **Real time clock.**
- **Internal Data Logger with computer Interface.**

📞 9891673394, 9716963092

✉ banditech07@gmail.com | info@banditechnology.com 🌐 www.indiamart.com/banditechnology