



Laboratory equipment > Homogenizers > Ov5 homogenizer





## OV5 Homogenizer

The ideal solution for dispersing, homogenizing and mixing biological tissues (cells, animal and vegetable tissues), pharmaceutical products, cosmetics and food products.

 CONTACT A SPECIALIST

€ REQUEST A QUOTE

 FIND A DEALER

share    




MODELS

ACCESSORIES

TECHNICAL DATA

DOWNLOAD

# Choose your Rotor/Stator





---

**STATOR:**

VS2C or VS3C coarse stator 20 or 30 mm

VS4C coarse stator 40 mm

VS2F or VS3F fine stator 20 or 30 mm

**ROTOR:**

VSR2 or VSR3 standard rotor 20 or 30 mm

VCR2 or VCR3 cutter rotor 20 or 30 mm

(for fibrous/stringy)

VMR2 or VMR3 mixing rotor 20 or 30 mm

(for solid/liquid)

VER2 or VER3 emulsifying rotor 20 or 30 mm

(for w/o or o/w emulsions)



Double clamp

cod. A00001301



Double Rod Mini H-Stand Set

cod. A00000383

**Stator rotor release kit**

cod. A00000279

**Strap clamp**

cod. A00000044

**VCR2 Cutter rotor 20 mm**

cod. A00000054

**VCR3 Cutter rotor 30 mm**

cod. A00000058



VER2 Emulsifying rotor 20 mm

cod. A00000056



VER3 Emulsifying rotor 30 mm

cod. A00000060



VMR2 Mixing rotor 20 mm

cod. A00000055



VMR3 Mixing rotor 30 mm

cod. A00000059





VS Shaft for homogenizer

cod. A00000047



VS2C Coarse stator 20 mm

cod. A00000048



VS2F Fine stator 20 mm


cod. A00000049



VS3C Coarse stator 30 mm

cod. A00000050



**VS3F Fine stator 30 mm**

cod. A00000051

**VS4C Coarse stator 40 mm**


cod. A00000052

**VSR2 Standard rotor 20 mm**

cod. A00000053

**VSR3 Standard rotor 30 mm**

cod. A00000057







VSS2CCR2 Dispersing tool

cod. A00000027



VSS2CMR2 Dispersing tool

cod. A00000028



VSS2CSR2 Dispersing tool

cod. A00000026



VSS2FCR2 Dispersing tool

cod. A00000031





cod. A00000023

cod. A00000032



VSS3CCR3 Dispersing tool

cod. A00000034




VSS3CMR2 Dispersing tool

cod. A00000036



VSS3CMR3 Dispersing tool



VSS3CSR3 Dispersing tool



VSS3FER3 Dispersing tool

cod. A00000037



VSS3FMR3 Dispersing tool

cod. A00000040



VSS3FSR3 Dispersing tool

cod. A00000038



VSS4CMR3 Dispersing tool

cod. A00000041



VSS5CSR4 Dispersing tool

cod. A00000046

Share



### VELP Scientifica Srl

Via Stazione, 16  
20865 Usmate (MB) - Italy  
Tel. + 39 039 628811

### VELP Scientific, Inc

155 Keyland Court, Bohemia,  
NY 11716 - U.S.  
Tel. +1 631 573 6002

### VELP China Co. Ltd.

Building 28, Lane 1333, Xinlong road  
Minhang district, Shanghai  
Tel. +86 18017557329



SCIENTIFICA

Copyright © 2020

**VELP Scientifica Srl**

Via Stazione 16 - 20865 - Usmate Velate (MB) - Italy

Tax Code: 06955700155 - VAT number: IT 00842180960

Company Registration Number MB: 06955700155 - REA number: MB-1129804 - Paid up  
share capital: € 500.000 fully paid.

[\[Privacy and Cookie policy\]](#) [\[Terms of Use\]](#)

